

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Account#: \_\_\_\_\_

Doctor: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State/Zip: \_\_\_\_\_ Telephone: \_\_\_\_\_

PATIENT NAME: \_\_\_\_\_

Dx: \_\_\_\_\_ Med Rec #: \_\_\_\_\_

Asymmetric Feet?  Yes  Male  
 No  Female Age: \_\_\_\_\_ Weight (required): \_\_\_\_\_

**For PROLAB  
office use only**

## CUSTOM ORTHOSES FOR SANDALS

ProLab manufactures custom orthoses for sandals with removable insoles (no thong-type).  
Sandals are not available through ProLab. **Sandals must be sent with all orders.**

### POLYPROPYLENE SHELL - All devices are direct milled (CD)

- Semirigid (**default**)  
or  
 3mm (1/8")  4mm (5/32")  5mm (3/16")  6mm (1/4")

### HEEL CUP DEPTH

- Match to shoe depth

### CAST FILL

- Minimum  R  L  
Standard  R  L  
Maximum  R  L

### MEDIAL HEEL SKIVE

- 2mm  R  L  
 4mm  R  L

### INVERT

- \_\_\_\_\_ ° Right (max 4°)  
\_\_\_\_\_ ° Left (max 4°)

### REARFOOT POST

- 0/0 motion poly post

### TOPCOVERS (to toes)

- Black**  Sport  Nylene  Leather/Poron  
**Gray**  Ultrasuede/Poron  
**Lt Brown**  Ultrasuede/Poron

### FOREFOOT EXTENSION

- Firm EVA beveled on device  
 Firm EVA Reverse Morton's Ext  R  L

### ADDITIONS

- Metatarsal Pad  R  L  
 Neuroma Pad  R  L \_\_\_\_\_ interspace

### SANDAL INFORMATION

**Sandals must be shipped with Rx**



**Sandal Type**

**Brand:** \_\_\_\_\_

**Color:** \_\_\_\_\_

**Shoe Size:** \_\_\_\_\_

### SPECIAL INSTRUCTIONS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

**DOCTOR'S  
SIGNATURE**  
(required) \_\_\_\_\_

### ORDER PROCESSING

(rush orders ship overnight)

- Rush 1 day in lab  
 Rush 3 day in lab

### SHIPPING

- Standard  
 FedEx Overnight  
 FedEx Ground  
 Ship to Patient  
(include address)

### REORDERS

**Do not fax - send sandals with Rx**

**Image#** \_\_\_\_\_  
(reference # from bottom of orthoses)

# CUSTOM ORTHOSES FOR SANDALS

## Descriptions, Policies, and Pricing

### ProLab Custom Orthoses for Sandals

ProLab manufactures custom functional orthoses for sandals with removable insoles (no thong-type). Orthoses for sandals are direct-milled polypropylene devices with a 0/0 rearfoot post. Orthoses match the width and heel cup depth of the sandal footbed. The topcover and forefoot extension matches the forefoot sandal depth (+/- 2mm). If you want the heel cup to be *higher* than the sandal depth, special instructions must be written. See *Custom Orthoses for Sandals Prescription Form* for available options (reverse side).

---

### Topcover Options

Topcovers are full length and provide a cushioning layer and a soft interface for comfort against bare skin.

#### Topcover Material Descriptions

<b>Sport (Black)</b>	Top layer synthetic suede, bottom layer high-rebound EVA
<b>Nylene (Black)</b>	Top layer nylon, bottom layer neoprene
<b>Leather/Poron (Black)</b>	Top layer leather, bottom layer Poron
<b>Ultrasuede/Poron (Gray or Lt Brown)</b>	Top layer Ultrasuede, bottom layer Poron

---

### Policies

- New Orders – send sandals, non-weightbearing casts, and Rx to ProLab
- Reorders – send sandals and Rx to ProLab
- Sandals must have a removable insole
- No thong-type sandals accepted
- Orthotic width and depth standard is within +/- 2mm of sandal footbed
- Due to the structure of sandals the orthoses may be visible. This is normal and there are no adjustments for visibility.

---

### Pricing

- Custom orthoses pricing applies plus a sandal fitting fee
- Shipping fees apply
- Rush orders available for additional cost: 1 Day in lab, 3 Days in lab
- No reorder discounts for sandal orthoses